



Spectradyne nCS1

Operation Manual



Contents

The nCS1 Instrument: Overview & Specifications.....	4
Specifications.....	4
Instrument Setup and Sample Prep.....	5
Setting up the nCS1	5
Cleaning cartridge.....	5
Control Sample and First Run.....	5
Loading the Cartridge	5
nCS1 Control and Data Acquisition Software	6
Priming the Cartridge.....	6
Data Acquisition	7
Turn on bias and sensing.....	7
Enter acquisition information.....	8
Acquire data (single acquisition mode).....	8
Assessing data quality.....	9
Acquire data (continuous acquisition mode):	9
The Auto-analysis Engine:.....	10
After Data Acquisition.....	10
Flush the cartridge:	10
Remove the cartridge:.....	11
nCS1 Shutdown Procedure	11
Preparing for long periods of inactivity	11
Drain fluid lines:	11
Remove the running buffer bottles and refrigerate	11
Turn off power to the instrument	11
Maintenance and Sample Preparation Guidelines.....	12
Buffer Maintenance	12
Cleaning Cartridge.....	12
Aborting a measurement run.....	12



Guidelines for Sample Preparation	12
Data Analysis and Viewer Software	14
Workflow Overview:	14
Data Viewer Software User Interface:	15
1. Select & Process Raw Data.....	15
Viewing and Selecting Raw Data for Analysis.....	16
Processing raw data.....	16
Viewing the Results of Raw Data Processing.....	17
2. Refine & Filter	18
2.1 Refine the diameter and concentration scaling of the CSD:	19
2.2 Filter the Set of Detected Particles to Exclude False-positives.....	20
3. Combine and Compare.....	22
3.1 Combine multiple Stats files into a single Combined file.....	22
3.2 Compare Multiple Analyses.....	22
More Information and Technical Support.....	24



The nCS1 Instrument: Overview & Specifications

The nCS1 is a powerful tool for obtaining detailed size and concentration information about particle distributions in liquids. The core technology utilized in the nCS1 is called Resistive Pulse Sensing (RPS), also referred to as the Coulter Principle, and is a truly orthogonal method to both light scattering methods such as dynamic light scattering or optical tracking as well as to resonant mass measurement techniques. A clean disposable cartridge targets a specific size range of particles for each analysis, eliminating concerns of contamination between samples.



Figure 1. The nCS1 instrument and cartridges.

Specifications

- nCS1 dimensions: 12" x 12" x 15" (30.5cm x 30.5cm x 38.1cm)
- nCS1 weight: 30 lbs (13.6 kg)
- Power requirements: North America: 110V, 5A, 60 Hz; Europe: 230V, 16 A, 50 Hz
- Nitrogen or CDA: 15 to 30 psi regulated supply
- Min. volume of analyte: 3 μ l
- Min. ionic strength of analyte: 1-10 mM
- Concentration range: $10^6 - 10^{12}$ particles per ml, cartridge dependent
- Running buffer: PBS buffer + 1% polysorbate 20 filtered to 20 nm is strongly recommended

A laptop computer is required to interface with the nCS1 through a standard USB link.



Instrument Setup and Sample Prep

Setting up the nCS1

Unpack: Remove all packaging from the instrument and position it on the benchtop.

Power supply: Use the power cable provided with the instrument to connect the nCS1 to a source of electrical power (North America: 110V, 5A, 60 Hz; Europe: 230V, 16 A, 50 Hz). The power input connection is located on the back panel. The fuse drawer also provides a means to set the instrument to run on 110V power (North America) or 230V (Europe). Align the arrows for the desired AC voltage on the fuse drawer and the power entry module.

Pressure supply: Connect the instrument to a source of clean inert gas (e.g., filtered nitrogen or air). The acceptable range of input pressures is 15-30 PSI. The pressure input is located on the back panel.

Load running buffer: Running buffer is used to fill auxiliary portions of the analysis cartridges during operation—it does not contact the sample to be analyzed until after measurement. Recommended running buffer is PBS + 1% Polysorbate 20 filtered at 0.02 um (see notes below). Add 25-50 mL to each of the two bottles labeled “Buffer” on the right side of the instrument.

USB connection: Use the USB cable provided with the instrument to connect the nCS1 to a computer with Spectradyn software installed. The USB port is located on the back panel.

Power on: Turn on the instrument with the switch located on the rear panel *before running the nCS1 software installed on the attached computer.*

Software: Run the Spectradyn nCS1 software on the computer. Note that if the instrument is powered off, the software will need to be restarted before next use.

Cleaning cartridge

After setting up the nCS1 for the first time, or after any significant period of inactivity, the cleaning cartridge routine should be performed using the reusable cleaning cartridge supplied with the instrument. Refer to page 5 of this guide for instructions.

Control Sample and First Run

A control sample is included with the nCS1. The control sample should be the first sample run after unpacking the instrument, and can be used at any time by users wishing to familiarize themselves with instrument operation. The control sample contains polystyrene nanoparticles in filtered PBS + 1% Tween 20. Guidelines for preparing your own sample for measurement can be found in the Maintenance and Sample Preparation Guidelines section on page 12 of this guide.

Loading the Cartridge

A video demonstrating cartridge loading is available on Spectradyn’s website: www.nanoparticleanalyzer.com.



- Pipette 3-5 microliters of the sample to be measured into the analyte reservoir on the analysis cartridge (**Figure 2**).
- Load the cartridge into the socket in the front of the nCS1. **IMPORTANT:**
 - ❖ The outer glass edge of the cartridge **MUST** be flush with the edge of the shelf on which it sits
 - ❖ The center indicator light (blue) and the cartridge engage button on the front panel of the instrument must be illuminated before proceeding.
- Ensure the nCS1 control software is running on the attached computer.
- Press the cartridge engage button on the front of the nCS1 to raise the cartridge into the instrument. The top indicator light (white) will illuminate to indicate that the cartridge is ready for priming.
- After this point the nCS1 and data acquisition process are controlled entirely through the Spectradyn Acquisition software on the computer.




Figure 2. Left panel shows the sample pipetted into the analysis cartridge. Right panel shows the front panel of the Spectradyn nCS1: On the left side is the cartridge insertion slot, and on the right are the cartridge engage button and indicator lights.

nCS1 Control and Data Acquisition Software

Priming the Cartridge

- Once the cartridge is loaded and engaged, it is ready for priming.
- Priming is an automated process in the nCS1—Press the “Prime” button in the acquisition software to begin the routine (see Figure 3).
- When the routine is complete, the software beeps and a message is displayed in the status window.
- After successful priming of the cartridge, the fill monitors should have approximately the pattern shown in Figure 3.
- The priming routine typically takes 1-2 minutes to complete. During the process, it is normal for fill monitors to vary between all levels and show different colors.




Figure 3. Fill monitors after successful priming of the cartridge.

Troubleshooting Tip:

- If the priming routine encounters an issue, one or more of the fill monitors will remain in a red color and the priming routine will not complete.
- After 5 minutes, a red fill monitor indicates that the corresponding portion of the cartridge was not successfully primed.
- If Fill Monitor 1 or 3 is red:
 - ❖ An issue with the running buffer is likely.
 - ❖ Check that sufficient running buffer remains in the reservoirs
 - ❖ Replace the running buffer and use the cleaning cartridge to flush the system with the new running buffer.
- If Fill Monitor 2 or 4 is red:
 - ❖ There may be an issue with the cartridge.
 - ❖ Try flushing the cartridge and restarting the Prime routine.
 - ❖ Try measuring the sample in a new cartridge.

Data Acquisition

Turn on bias and sensing

- Refer to Figure 4
- The -V and +V bias entries will be preset to -2V and +2V respectively.
- Click the Bias ON/SET button to turn on the bias




Figure 4. Bias and sensing frame. The “Bias ON/SET” button is used to turn on the bias and/or set the bias to the input values.

- After turning on the bias, the Sensor DC level should remain approximately centered on the scale.
- During normal operation the Sensor DC level should be in the range of -0.400 V – +0.400 V, and the color of the scale should be green.
- If the scale moves outside the normal range, run the “Clear Constriction” routine described below.

Enter acquisition information

Description of fields (* denotes required, see **Figure 5**):

- **Mold ID***: Enter the middle section of the cartridge box label in this field, e.g., P1.2.3, or E16.1.2.
- **CalPs**: Enter the diameter (in nm) of control particles mixed in with the sample. For example, for a mix of 94 nm and 150 nm control particles, enter, “94, 150” in the CalPs field. If the sample contains no control particles, the entry may be left blank.
- **CalP Conc**: Enter the concentration of calibration particles (in particles/mL). For a mix of 94 nm and 150 nm control particles, each at 3E9/mL, enter “3E9,3E9” in the CalP Conc field. This entry may be left blank.
- **Sample**: Enter a sample description, this will be stored with the data file.
- **Comment**: Enter any comments, these will be stored with the data file.




Figure 5. Data acquisition frame. A valid Mold ID is required before acquisition functions are enabled. Use the “Acq+Plt+Save” button to acquire, plot, and save data directly to file.

Acquire data (single acquisition mode)

- Set the acquisition time (in seconds) using the Acq. (s) dropdown list.
- Click “Acq+Plt+Save” to Acquire, Plot, and Save the raw data with a single click.
- Click “A&P” to acquire and plot the raw data without saving
- Subsequently clicking the “Save” button will save the last acquisition.



Data Saving Tip:

- Raw data files are saved in a dated subdirectory of the Spectradyne data directory, C:\Spectradyne\Data\. For example, if data is acquired on December 1 2015, a file will be saved in the directory: C:\Spectradyne\Data\2015_12_01 by default.
 - Files are named with a date and a sequential number starting at 0 that is incremented with each save. For example, a raw data file acquired on December 1 2015 might be named 20151201_003.h5.
- The current acquisition can be interrupted using the “Stop Acq.” Button. Partial acquisitions are not saved, and no new data will be displayed.

Assessing data quality




Figure 6. Exemplary raw data: Smooth baseline with particle detection events appearing as downward peaks.

- Raw data should have the form shown in **Figure 6**—it has a smooth baseline with particle detection events appearing as downward spikes from the baseline.
- Drifting of the baseline is acceptable.
- In some situations, there may be occasional steps in the baseline. These steps may be acceptable, depending on their severity.
- If the raw data shows too many steps, or poor display of control particles, click on the “Clear Constriction” button (Figure 3), wait for completion, then reattempt acquisition.
- In single acquisition mode, the user may flag the most recent acquisition as “Good” by pressing the “ON” button in the Last Acquisition frame. Flagged acquisitions are highlighted for rapid identification in the Viewer software (refer to the Data Analysis and Viewer Software section on page 13 of this guide).

Acquire data (continuous acquisition mode):

- Under normal operation all acquisitions are high quality, and the user may enter continuous acquisition mode by pressing the “Continual Acq.” button.



- In this mode, acquisitions of length specified by the “Acq. Time” dropdown, are made repeatedly until the user presses the “Stop Acq.” button.
- Acquisitions made in continuous mode are flagged as “Good” by default; the user may un-flag the latest acquisition by pressing the “OFF” button in the Last Acquisition frame if necessary.

The Auto-analysis Engine:

- Flagged acquisitions are placed in a queue for automatic processing into Spectradyne Stats files by the Auto-analysis Engine.
- The time required for the auto-analysis engine to process raw data depends on the concentration of particles in the sample, but is typically near real-time.
- As a result, the measured particle size distributions can be examined in parallel with the acquisition process using the Viewer software, and used to evaluate the quality and quantity of data acquired to that point.
- The auto-analysis engine uses standard default parameters for processing (refer to the section, “Processing raw data” on page 16 of this guide for details):
 - Peak Detection Threshold: 3.0 x the standard deviation of the baseline noise.
 - Fit Width: 4.0 x the mean transit time of high signal-to-noise particle events.

How much data is enough?

- The 1- σ error reported for a concentration measurement in a given size range is derived from statistical error in counting particle detection events.
- Particle detection events are Poisson-distributed in time; therefore, the error scales as $1/N^{1/2}$ (with small corrections for $N < 5$).
- Example: $N = 100$ particles must be counted in a given size range to obtain a concentration measurement on that range having standard error of 10 %.
- The particle concentration (with error) measured over any user-defined size range is easily obtained from the particle size distribution using the tools on the “Quantification” tab of the Viewer software.

After Data Acquisition

Flush the cartridge:

- After sufficient data has been acquired, flush the cartridge.
- Flushing is an automated process in the nCS1—Press the “Flush” button in the acquisition software to begin the routine (see Figure 3)
- When the flush routine is complete, the software beeps and a message is displayed in the status window.
- After successful flushing of the cartridge, the fill monitors should all be close to zero.
- Completion of the flush routine typically requires 1-2 minutes. During the process, it is normal for fill monitors to vary between all levels and show different colors.



Remove the cartridge:

- After completion of the flush routine, lower the stage to remove the cartridge by pressing the button on the front panel of the instrument.
- A new cartridge may now be loaded for analysis.

nCS1 Shutdown Procedure

Preparing for long periods of inactivity

Before long periods of inactivity (e.g., longer than about 1 week), the cleaning cartridge routine should be performed using **deionized water in place of the normal running buffer**.

- Remove the buffer bottles from the instrument:
 - Drain fluid lines using the procedure described below.
 - Detach the buffer bottles from the fluid lines using the quick-connect tube fittings and take the bottles out of the instrument.
- Remove the running buffer from the bottles and refrigerate if intended for future use.
- Add approximately 25-50 mL of deionized water to each of the two buffer bottles, and replace in the instrument.
- Run the cleaning cartridge routine as described on page 5 with DI in place of running buffer.
- Wait for the status message to indicate that the process is complete, then continue with the shutdown instructions below, starting by draining the fluid lines.

Drain fluid lines:

- In the nCS1 control software, click the “More” button to reveal the Utilities frame.
- Click the “Instrument Shutdown” button to drain the fluid lines and prepare the instrument for shutdown.
- Wait for the status message to indicate that the process is complete.

Remove the running buffer bottles and refrigerate

- The running buffer bottles should be refrigerated during periods of inactivity (e.g., overnight) to minimize algae growth.
- Detach the buffer bottles from the supply lines using the quick-connect tube fittings and store in the refrigerator until next use.
- Running buffer should be replaced with fresh filtered buffer approximately every two weeks

Turn off power to the instrument

- It is recommended to power off the instrument during long periods of inactivity (e.g., overnight). Because the nCS1 detects particles electrically (not optically), the instrument requires no warm up time.
- Use the power switch at the back of the instrument to power down.
- The power indicator light (green) on the front of the instrument may take a few seconds to turn off.



Maintenance and Sample Preparation Guidelines

Buffer Maintenance

Refrigerate the running buffer bottles

- The running buffer should be refrigerated during long periods of inactivity (e.g., overnight) to minimize algae growth.
- Running buffer should be replaced with fresh filtered buffer approximately every two weeks.

Cleaning Cartridge

After long periods of inactivity (e.g., 1-2 weeks) the cleaning cartridge should be used to flush the system with clean running buffer before analyzing a new sample.

- Load 25-50 mL filtered running buffer into each buffer bottle on the right side of the instrument (see Setting up the nCS1, page 5).
- Insert the cleaning cartridge provided with the nCS1 and raise it into the instrument
- In the nCS1 control software, click the “More” button to reveal the Utilities frame.
- Click the “Run Cleaning Cartridge” button to begin the automated cleaning process.
- Wait for the process to complete, then remove the cleaning cartridge
- The nCS1 is ready for normal operation using a new analysis cartridge.

Aborting a measurement run

The best way to finish using a cartridge is to complete the automated flush routine as described in the section entitled, “After Data Acquisition” on page 10. In rare cases, if the flush routine fails, it may be necessary to eject the cartridge manually and reset the instrument. Follow these steps:

- Try the flush routine first. If it fails to complete then proceed to the next steps.
- Cancel the flush routine using the “Cancel” button (see Figure 3).
- Click the “More” button to reveal the Utilities frame.
- Click the “Abort Cartridge” button and follow the instructions provided.
** Note: the cleaning cartridge routine (see page 12) should always be run after aborting a measurement run.

Guidelines for Sample Preparation

Some optimization of sample preparation may be required when analyzing a sample for the first time. Once the best sample prep for the analyte in question is determined, that procedure can be used for future analysis of that analyte and its variations. The parameters below are suggested targets for first-run samples.

Optimal sample parameters:

Concentration:	10^9 particles/ml below 300 nm, 10^8 particles/mL above 300 nm
Control beads:	Add control beads to final concentration $\sim 10^9$ /ml below 300 nm, $\sim 10^8$ /ml above 300 nm. Recommended size $\sim 30\%$ below the cartridge max size.
Ionic strength:	$\sim 1\text{S/m}$ (e.g., $\sim 100\text{ mM NaCl}$)



Surfactant: 1% polysorbate 20 or equivalent

Notes:

- Sample diluent and control particles:
 - ❖ If ionic strength or surfactant in the diluent is a concern, then lower values can be used.
 - ❖ For first-run analyses, Spectradyne highly recommends mixing control particles of a known size (typically NIST-traceable polystyrene beads) with the analyte.
 - ❖ Keep in mind that the sample diluent may be a source of nanoparticles!
- Running buffer:
 - ❖ **Important:** Running buffer does not contact the analyte sample until after its particle distribution has been analyzed. The user should not be concerned that ionic strength or surfactant content of the running buffer will affect the measurements.
 - ❖ Running buffer should be refrigerated during periods of inactivity (e.g., overnight)
 - ❖ Running buffer should be replaced approximately every 2 weeks in order to avoid contamination of the fluid lines from algae, etc.



Data Analysis and Viewer Software

Workflow Overview:




Figure 7. Overview of data processing workflow.

Three workflow steps for data analysis

1. Select & Process:
 - Raw data acquisitions are selected and processed into Spectradyne Stats files. Likely this step will have been completed by the Auto-analysis Engine at the time of acquisition (described on page 7), however a manual process is available for customizing the choice of raw data and processing parameters.
 - This step yields a quantitative Concentration Spectral Density (CSD™) distribution.
2. Refine & Filter:
 - Any necessary refinement to the diameter scaling is performed (e.g., to account for possible buffer effects)
 - Low-confidence particle measurements are eliminated by data filtering (e.g., false positives near the small particle detection limit)
 - This step yields the final CSD for the collection of measurements.
3. Combine & Compare:
 - The collection of refined Spectradyne Stats files is saved as a single combined file.
 - The CSDs of different samples are quantitatively compared at the combined file level.



Data Viewer Software User Interface:




Figure 8. Overview of the Spectradyne Viewer software user interface.

1. Select & Process Raw Data

In a typical workflow, data to be used in analysis will have been flagged “Good” at the time of acquisition, and processed into Spectradyne Stats files by the Auto-analysis Engine (see page 10). In such a scenario the user may skip to the section entitled, “Viewing the Results of Raw Data Processing” on page 17, below. If the raw data has not already been processed, or if the user wishes to employ non-standard processing parameters (i.e., a user-defined peak fit width or peak threshold), then the manual method for processing raw data described next should be used.



Viewing and Selecting Raw Data for Analysis

Navigate to data location

- Open the Data Viewer software.
- In the Select/Quickview frame (**Figure 8**), use the folder icon to the right of the current path to navigate to the location of the raw data.

Inspect raw data

- Check one of the plot options in the Raw Data Plot Parameters frame:
 - ❖ “Plot Raw” – show the data as it was collected
 - ❖ “Plot Polished” – show signal-processed data in which particle events will be measured.
- Select a file from current folder in the Select/Quickview list box—the raw data will be displayed.
- The “Prev” and “Next” buttons frame navigate through the files in the current directory for rapid inspection of raw data.

Note: Select/Quickview Frame

- Files flagged as “Good” are highlighted with blue text.
- The flag can be toggled for selected files using the “Marked as Good Data” checkbox.
- The Select/Quickview list can be filtered to show subsets of the folder contents using the checkboxes below the list (e.g., raw files only).
- Metadata for any selected file is shown in the Active File Info frame.

Load raw data files for processing

- After identifying the raw data to be analyzed, load the files into the Currently Loaded file list at right using one of two methods:
 1. Use the load button (➡) to load selected files from the Select/Quickview list.
 2. Use the remote load button (📁➡) to select files from another folder with a dialog window.
- The Currently Loaded file list may be modified with the control buttons at the bottom of the frame.

Processing raw data

Each raw data file is analyzed to count and measure particles detected during acquisition. The result of this analysis is a Spectradyne Stats file that contains detailed information about all particles detected during the acquisition as well as the acquisition itself.

Procedure for processing raw data

- The Process Raw Data tab of the Data Viewer contains all the controls for processing raw data, as well as relevant measured values for the selected raw data file.
- The default processing parameter values are appropriate for use in most cases. These default values are the same as those used by the Auto-analysis Engine during acquisition:
 - ❖ Peak Detection Threshold: 3.0 x the standard deviation of the baseline noise.
 - ❖ Fit Width: 4.0 x the mean transit time of high signal-to-noise particle events.
- The above parameters can be customized for more in-depth analysis using the entries in the Processing Parameters frame.
- To analyze the raw data, use the controls in the Processing Actions frame:



- ❖ “Process Selected”: Processes only the currently selected raw data file from the Currently Loaded list.
- ❖ “Process All”: Processes all raw files in the Currently Loaded files list using the same (adaptive) parameters.
- ❖ “Stop”: Interrupts processing of the current raw file with no partial data saved and stops batch processing.

Note: Stats Files

- A Stats file shares the same base filename as the raw file from which it was generated, with an additional extension of the form `_ssNN`.
 - In the extension, NN is an integer that increments each time the raw file is processed.
 - Stats files are saved in the same folder location as the raw file from which they are generated.
 - For full functionality, a Stats file should be kept in the same folder as its parent file.
- Processed raw files are replaced by the newly generated Stats file in the Currently Loaded list.

Viewing the Results of Raw Data Processing

Ensure that the Stats files of interest are loaded into the Currently Loaded files list.

Concentration Spectral Density (CSD):

- The nCS1 is capable of very high resolution size and concentration measurements of particles in a liquid, and is well suited for characterizing and comparing subpopulations of particles in polydisperse mixtures.
- A useful quantity is therefore the concentration of particles in the sample over a given size range, so that different size populations in a mixture can be compared directly.
- The Concentration Spectral Density (CSD) is a measure of particle concentration (particles/mL) per unit diameter (nm). The CSD is analogous to a probability density distribution.
- The integral of a CSD over a given size range then yields the concentration of particles (in units of particles/mL) in that range.

Aggregate Data View: Generating the CSD

- Select a plot mode from the QuickView Processed Data frame:
 - ❖ “Single” – Display results from only the selected Stats file in the Currently Loaded file list.
 - ❖ “Multi” – Display results from all Stats files in the Currently Loaded file list separately on the same plot.
 - ❖ “Combined” – Combine results from all Stats files in the Currently Loaded file list as if they derived from a single long acquisition.
- Plot the CSD using one of two available “CSD” buttons:
 - ❖ In the QuickView Processed Data frame (always visible)
 - ❖ In the Aggregate Views frame of the Data Visualization tab.



Single-particle Views: Particle events in raw data

- Navigate to the Data Visualization tab in the DataViewer.
- Use the “Peaks in Raw Data” button to show all peaks detected in the stats file that is selected in the Currently Loaded list.
 - ❖ Use the tools in the plot window to zoom in on individual peaks.
 - ❖ Clicking on the peak (at sufficiently high zoom) displays the details of the peak measurement (**Figure 9**).
 - ❖ Ctrl + Click selection allows concurrent display of multiple peak measurements.




Figure 9. Inspecting individual peaks detected in the raw data.

Single-particle Views: Parametric Scatter Plots

- The “Scatter” button in the Data Visualization tab generates a scatter plot of particle diameter vs the parameter selected in the dropdown list (e.g., transit time) for each detected particle event.
- The plot mode selected in the QuickView Processed Data frame determines whether data from a single Stats file or multiple Stats files will be displayed.

Note: Refine Data Group

Any outlier Stats files should be eliminated from the loaded file list at this point in the analysis workflow.

- Generate a CSD plot in multi-mode to show the CSDs for all Stats files in the Currently Loaded list.
- CSDs corresponding to measurements of the same sample should be consistent with each other—any outliers should be removed from the Currently Loaded list using the red “X” button below the list.

2. Refine & Filter

At this point in the data analysis workflow, multiple measurements of the same sample will likely have been processed to generate Stats files in the Currently Loaded files list. These Stats files will contain measurements of many individual particles in the same sample under the same measurement conditions, and can therefore be considered as a single analysis of the sample for downstream processing. The



Spectradyne Viewer provides efficient tools for working with multiple Stats files as a group, and after the Refine & Filter workflow step the multiple files are combined into a single entity for comparison with other sample measurements.

Note: Treating multiple datasets as one for Refine and Filter

- In the QuickView Processed Data frame of the Spectradyne Viewer, use the Combined plot mode to display data from all the loaded Stats files as a single dataset.

The two main goals of this step in the workflow are:

1. Refine diameter and concentration scaling as necessary.
2. Filter the particle set to keep only high-confidence detection events and eliminate false-positives.

2.1 Refine the diameter and concentration scaling of the CSD:

Particles are detected and measured in an analysis cartridge that is cast from a master mold. Different cartridges may be cast from different molds, and the analysis software automatically compensates for variation between cartridges using a mold-specific calibration factor. This factor establishes the correct scaling for the diameter axis of the CSD in standard measurement conditions. Certain factors, (e.g., composition of the sample diluent), may cause slight deviations from the default calibration, and the user may wish to refine the diameter scaling to align known features of the CSD (e.g., a peak of control particles) with the appropriate diameter or concentration.

Procedure for refining scaling

- Select the Quantification tab in the Spectradyne Viewer
- In the Refine Scaling frame, use the Manual Scaling Refine entries to scale the diameter and concentration. The value of the entry multiplies the CSD; for example, a Diameter refine factor value of 1.02 would ‘stretch’ the diameter scale by a factor of 1.02, and a population of particles originally centered at 100 nm would shift to a new center of 102 nm.
- Occasionally, more challenging samples may require adjustment of the concentration scaling in an analogous way.

Note: Refining the diameter scaling using a known feature in the CSD

The diameter scaling of a CSD can quickly be refined using a known feature in the CSD (e.g., a population of calibration particles added to the sample that exhibit a clear peak in the particle size distribution):

- Fit a peak in the CSD to a Gaussian curve using the steps described in the section, “Quantifying Particle Size Distributions” on page 22 of this guide.
- The “Rescale Last Fit To:” button will automatically refine the diameter scaling so that the mean of the Gaussian fit coincides with the diameter specified by the user in the entry adjacent to the button.



2.2 Filter the Set of Detected Particles to Exclude False-positives

Particle detection events manifest as peaks in the raw data. The algorithm used to detect and characterize these peaks employs a thresholding method: Any excursion of the normalized output signal crossing above the peak detection threshold is considered to be a possible peak and measured (see **Figure 10**). The metrics for each peak found in this way include quality of fit parameters (i.e., signal-to-noise), as well as physical parameters of the transit event such as transit time through the sensing constriction and symmetry of the peak shape.

As the peak detection threshold approaches the baseline noise, some false positives (i.e., peaks in the raw data that do not likely represent real particles) naturally are included in the Stats file data set (see **Figure 10**). The goal of this step in the workflow is to eliminate false positives from the data by excluding particle detection events whose measured peak parameters fall outside the normal range for true particles. The signal-to-noise, peak symmetry, and transit time parameters are all available for use as filter parameters.




Figure 10. The peak detection threshold indicated by the dashed line is close enough to the baseline noise that two false positives are detected near 200 us and 450 us in the local time window; a real particle event is clearly detected at 750 us. A red point highlights the peak value, and the red line shows the fit for each peak. Peak-filtering tools are available in the Viewer for identifying and removing false-positives in the data.

Identifying False-Positives Among the Measured Peaks

False-positives generally are characterized as follows (example shown in **Figure 11**):

- Diameter:
 - Diameter near or below the small size limit of detection of the cartridge.
- Transit time:
 - A group of small diameter particle events having a broad distribution of transit times.
- Symmetry:
 - Symmetry outside the range of 0.2 – 4.0 (approx.)
- Signal-to-noise (S/N):
 - $S/N < 15$ (cartridge dependent).
 - A general sense for an appropriate S/N cutoff can also be obtained by inspecting peaks in the raw data (**Figure 10** and “Viewing the Results of Raw Data Processing,” page 17).




Figure 11. Representative Scatter plot of Diameter vs. Transit Time showing two distinct populations of true particle detection events (red) and false-positives excluded from the data set (blue). False positive detection events are characterized in part by having diameter close to the minimum detectable size for the cartridge (here approx. 130 nm) and a broad distribution of transit times.

Procedure for removing false-positives from the set of detected peaks

- Load the desired Stats files into the Currently Loaded files list
- Navigate to the Peak Filtering tab of the Spectradyne Viewer
- Use the appropriate action buttons in the Filters frame of the Peak Filtering tab to apply and remove existing filters as desired from files in the Currently Loaded file list.
- To create a new filter, use the “+” button at the bottom of the available filters list
- In the New Filter dialog box:
 - Choose one of two options for type of filter:
 - Linear – Filter bounds have constant value of the non-diameter parameter and include an infinite range of diameters
 - Polygon – Filter has 2D, user-defined polygonal boundaries similar to a gate commonly used in analysis of flow cytometry results.
 - Select the filter parameter of choice from the dropdown (e.g., Transit Time)
 - Specify a descriptive name for the filter
 - Specify the limits of the filter in one of two ways:
 - Enter min and max values for the non-diameter parameter in the entries provided (linear filters only)
 - Specify the bounds of the filter directly on a scatter plot using the mouse
 - Click “Create” to add the filter to the list of available filters.
- The effects of peak filtering will be reflected in any of the viewing methods available (e.g., CSD plots, parametric scatter plots). In Single and Combined viewing mode, both the raw and filtered datasets are displayed in the scatter plots.



The refined and filtered Stats files can now be saved individually or combined into a single Combined file format as described in workflow step 3, below. To save one or more Stats files with the applied filters and scaling, navigate to the Data Output tab of the Viewer and use the controls in the Save frame.

- Saving a Stats file creates a new file with an incremented extension number, and does not overwrite existing Stats files present in the current working folder.

3. Combine and Compare

At this point in the workflow, the Currently Loaded file list typically contains multiple Stats files that together represent a polished analysis of the particle content of a single physical sample. The goals of this workflow step are to:

1. Combine the multiple Stats files into a single entity for simpler manipulation.
2. Use the Viewer's quantitative analysis tools to compare the measurement to others.

3.1 Combine multiple Stats files into a single Combined file

Combine Procedure

- Load the desired Stats files into the Currently Loaded file list.
- Navigate to the Data Output tab of the Spectradyne Viewer.
- In the Combine Stats frame, use the "All Loaded" button to combine all loaded Stats files into a single Combined file.
- A new file will be generated in the current working directory with file name of the form, "YYYYMMDD_ccNN.h5", where NN is an incrementing number unique for each day.

3.2 Compare Multiple Analyses

- Load the analyses of interest into the Currently Loaded file list.
- In the QuickView Processed Data frame of the Spectradyne Viewer, select Plot Mode Multi.
- Use the tools in the Data Visualization and Quantification tabs of the Viewer to make quantitative comparisons of the analyses.
- Scatter and CSD plots will display separate results from each loaded analysis on the same axes, and quantitation tools will report results for each analysis separately.

Quantifying Particle Size Distributions – Concentration in the CSD

- The concentration of particles over a certain size range is readily obtained after plotting the CSD (**Figure 12**):
 - ❖ In the Quantification tab of the Viewer, click the "Select Range" button.
 - ❖ Click and drag the mouse cursor across the CSD over the size range of interest.
 - ❖ Click "Measure Concentration" to report the concentration over the selected size range.
 - ❖ The concentration measurement is displayed on the plot.

Quantifying Particle Size Distributions – Fitting Gaussian Peaks in the CSD

- Peaks in the CSD can be fit to a normal distribution to obtain the mean and width of the peak (**Figure 12**):
 - ❖ In the Quantification tab of the Viewer, click the "Select Range" button.



- ❖ Click and drag the mouse cursor across the CSD over the size range of interest.
- ❖ Click “Fit to Gaussian” to fit the selected region to a normal distribution.
- ❖ Fit results and the concentration of particles over the fit range are displayed on the plot

The above quantification methods perform the same functions when plotting multiple CSDs (multi-mode). The CSD corresponding to the file selected in the Currently Loaded Files dropdown is the one that is measured.




Figure 12. Quantification of a CSD plot. The peak is fit to a Gaussian distribution and the concentration of particles in the range of 200-300 nm is measured directly.

Note: Exporting Stats data to other file formats

Navigate to the Export frame in the Data Output tab of the Viewer.

- Two formats are available for export:
 - Microsoft Excel format (.xlsx)
 - MS DOS CSV format (.csv)
- When exporting to Excel format (“Excel” button), the file selected in the Currently Loaded files list is exported to two sheets in a single file named <original stats filename>.xlsx:
 - Sheet 1 contains the metadata associated with the file (e.g., acquisition parameters, processing parameters, etc.).
 - Sheet 2 contains the CSD data and a chart showing the CSD.
- When exporting to CSV format (“CSV” button), the selected file is exported to a single sheet in a file named <original stats filename>.csv. CSD data is placed ‘below’ the metadata in the sheet.
- The “All...” button saves all the stats files in the Currently Loaded files list to a single file.
 - A dialog window permits the user to choose the export format and the filename.



More Information and Technical Support

More detailed information on the operation of the nCS1 and its software can be obtained by sending email to support@spectradynellc.com.

At Spectradyne we are eager to help with your particle analysis. Please contact us for help with specific measurement challenges and for help with the analysis of specific particle populations.



info@SpectradyneLLC.com
www.SpectradyneLLC.com